



CODAN Syringe pumps

A616S TCI

A616S InCare

List of qualified syringes



List of qualified syringes

A616S InCare / A616S TCI

Brand	Name	Nominal Volume
Alaris		
Alaris	IVAC	50 ml
BD		
BD	Plastipak	10 ml
BD	Plastipak	20 ml
BD	Plastipak	30 ml
BD	Plastipak	50 ml
BD ¹	Perfusion	50 ml
B. Braun		
B. Braun	Omnifix	10 ml
B. Braun	Omnifix	20 ml
B. Braun	Omnifix	30 ml
B. Braun	Omnifix	50 ml
B. Braun	Perfusor	20 ml
B. Braun	Perfusor	50 ml
B. Braun	Inkjet	20 ml
CODAN		
CODAN		10 ml
CODAN		20 ml
CODAN		30 ml
CODAN		50 ml
Fresenius		
Fresenius	Injectomat	50 ml
Fresenius	Heparin	30 ml
Fresenius-Kabi	Kabifill PFS 50	50 ml
Fresenius-Kabi	Kabifill PFS 20	20 ml
Once		
Once		10 ml
Once		20 ml
Once		30 ml
Once		50 ml
Penta		
Penta		10 ml
Penta		20 ml
Penta		30 ml
Penta ²		60 ml
Terumo		
Terumo		10 ml
Terumo ³		20 ml
Terumo		30 ml
Terumo ³		50 ml

¹ Flow rate accuracy +/- 2% for rate \geq 5ml/h.

² Do only support infusion at max 55 ml.

³ Flow rate accuracy +/- 2% for rate \geq 5ml/h.



Manufacturer

CODAN ARGUS AG

Oberneuhofstrasse 10 · 6340 Baar, Switzerland

Tel. +41 (0) 41 · 785 09 44 · Fax. +41 (0) 41 · 785 09 40

www.codanargus.com

The decisive connection

